

Product Description SALSA[®] MLPA[®] Probemix P357-A3 MODY Mix 2

To be used with the MLPA General Protocol.

Version A3

For complete product history see page 10.

Catalogue numbers:

- P357-025R: SALSA MLPA Probemix P357 MODY Mix 2, 25 reactions.
- **P357-050R:** SALSA MLPA Probemix P357 MODY Mix 2, 50 reactions.
- P357-100R: SALSA MLPA Probemix P357 MODY Mix 2, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see www.mrcholland.com).

Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mrcholland.com.

Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mrcholland.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

General information

The SALSA MLPA Probemix P357 MODY Mix 2 is a **research use only (RUO)** assay for the detection of deletions or duplications in the *KLF11*, *NEUROD1*, *PAX4*, *CEL*, *INS*, *PDX1* and *HNF1B* genes, which are associated with Maturity-Onset Diabetes of the Young.

Maturity-Onset Diabetes of the Young (MODY) is a distinct form of non-insulin-dependent diabetes mellitus (NDDM), also known as type II diabetes. MODY has an autosomal dominant inheritance and it generally develops in individuals under 25 years of age. Approximately 5% of all diabetes patients suffer from MODY. As described in more detail below, 14 forms of MODY have now been identified with the highest prevalence for MODY 2 and 3. Each form of MODY has been associated with one gene.

This P357-A3 MODY Mix 2 probemix contains probes for the *PDX1*, *HNF1B*, *NEUROD1*, *KLF11*, *CEL*, *PAX4* and *INS* genes. This probemix is therefore specific for MODY 4-10. The widely used probemix P241 MODY Mix 1 contains probes for the *HNF4A*, *GCK*, *HNF1A* and *HNF1B* genes involved in MODY types 1-3 and 5, the more frequent causes of MODY. For these genes a clear relation has been described between copy number changes of the genes and the MODY phenotype or the renal cysts and diabetes (RCAD) syndrome.

The MODY genes

- **MODY 1** is a result of defects in the Hepatocyte nuclear factor-4-alpha gene (*HNF4A*) on chromosome 20q13. This gene comprises 12 exons and spans about 30 kb of genomic DNA. The HNF4A protein regulates the expression of *HNF1A*. MLPA probes for the *HNF4A* gene are present in the P241 probemix.
- **MODY 2** is caused by mutations in Glucokinase gene (*GCK*) on chromosome 7p13. This gene comprises 12 exons and spans about 45 kb of genomic DNA. MLPA probes for *GCK* are present in the P241 probemix.
- **MODY 3** is caused by defects in the HNF1 Homeobox A gene (*HNF1A*) which comprises 10 exons and spans 24 kb of genomic DNA on chromosome 12q24. MLPA probes for *HNF1A* are present in the P241 probemix.
- MODY 4 has been linked to defects in the Pancreas/Duodenum homeobox protein 1 gene (*PDX1*). The *PDX1* gene (two exons) spans ~6 kb of genomic DNA and is located on chromosome 13q12.2, ~27 Mb from the p-telomere.

- MODY 5 have been associated with the HNF1 Homeobox B gene (HNF1B). The HNF1B gene (nine exons) spans ~59 kb of genomic DNA and is located on chromosome 17q12, ~33 Mb from the p-telomere. The HNF1B probes in this P357 probemix detect the same sequences as the HNF1B probes in P241.
- MODY 6 has been linked to defects in the neuronal differentiation 1 gene (*NEUROD1*). The *NEUROD1* gene (two exons) spans ~4 kb of genomic DNA and is located on chromosome 2q31.3, ~182 Mb from the p-telomere.
- **MODY 7** is caused by mutations in the Krüppel-Like Factor 11 gene (*KLF11*) on chromosome 2p25.1. The *KLF11* gene (four exons) spans ~11 kb of genomic DNA and is located on chromosome 2p25.1, ~10 Mb from the p-telomere.
- MODY 8 has been associated with defects in Carboxyl-ester lipase gene (CEL). The CEL gene (11 exons) spans ~10 kb of genomic DNA and is located on chromosome 9q34.2, ~135 Mb from the p-telomere. Due to sequence similarity with the CEL pseudogene CELP, only probes for exons 2 to 7 of CEL are included. Compared to CEL, CELP lacks a 4.8-kb fragment containing exons 2 to 7. The CEL gene is located in a CNV region (http://dgv.tcag.ca/dgv/app/home). We observed duplications of this gene in DNA samples of several healthy individuals.
- **MODY 9** is caused by defects in the Paired box 4 gene (**PAX4**). The PAX4 gene (12 exons) spans ~8 kb of genomic DNA and is located on chromosome 7q32.1, ~127 Mb from the p-telomere.
- MODY 10 has been linked to mutations in the insulin gene (*INS*). The *INS* gene (three exons) spans ~1.4 kb of genomic DNA and is located on chromosome 11p15.5, ~138 Mb from the p-telomere. Please note that *INS* gene mutations that cause MODY are all heterozygous missense mutations with a dominant-negative mode of action; a heterozygous deletion of the *INS* gene might not result in diabetes. The ME033 probemix contains *INS* probes that target the same sequences as the *INS* probes in this P357 probemix.
- **MODY 11** has been linked to defects in the B lymphoid tyrosine kinase gene (*BLK*). *BLK* has 13 exons and spans about 71 kb of genomic DNA on chromosome 8p23.1. No probes for this gene are included in P241 or P357.
- **MODY 12** has been linked to defects in the ATP binding cassette subfamily C member 8 gene (*ABCC8*). The *ABCC8* gene consists of 39 exons and spans 84 kb of genomic DNA on chromosome 11p15.1. Probes for this gene are included in the P117 probemix.
- **MODY 13** has been linked to defects in the potassium inwardly rectifying channel subfamily J member 11 gene (*KCNJ11*). The *KCNJ11* gene consists of one exon and spans 3.4 kb of genomic DNA on chromosome 11p15.1. Probes for this gene are included in the ME033 probemix.
- **MODY 14** has been linked to defects in the gene encoding adaptor protein, phosphotyrosine interacting with PH domain and leucine zipper 1 (*APPL1*). The *APPL1* gene consists of 22 exons and spans 46 kb of genomic DNA on chromosome 3p14.3. No probes for this gene are included in P241 or P357.

More information is available at https://www.ncbi.nlm.nih.gov/books/NBK500456/.

This SALSA MLPA probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.

Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene For NM_ mRNA reference sequences: http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide Locus Reference Genomic (LRG) database: http://www.lrg-sequence.org/

Exon numbering

The *KLF11*, *NEUROD1*, *PAX4*, *CEL*, *INS*, *PDX1* and *HNF1B* exon numbering used in this P357-A3 MODY Mix 2 product description is the exon numbering from the NG_017199.1, LRG_1119, NG_012848.1, NG_016394.1, NG_007114.1, NG_008183.1 and NG_013019.2 sequences, respectively. The *PAX4* exon numbering has changed; the exon numbering used in previous versions of this product description can be found in between brackets in Table 2c. The exon numbering of the NM_ sequence that was used for determining a probe's ligation site does not always correspond to the exon numbering obtained from the NG or LRG sequences. As changes to the databases can occur after release of this product description, the NM_ sequence and exon numbering may not be up-to-date.

Probemix content

The SALSA MLPA Probemix P357-A3 MODY Mix 2 contains 45 MLPA probes with amplification products between 131 and 490 nucleotides (nt). This includes four probes for the *KLF11* gene (one for each exon), two probes for the *NEUROD1* gene (one for each exon), nine probes for the *PAX4* gene (for exons 4 to 12), five probes for the *CEL* gene, three probes for the *INS* gene (one for each exon), two probes for the *PDX1* gene (one for each exon) and ten probes for the *HNF1B* gene (one for each exon and two for exon 4). In addition, ten reference probes are included that detect autosomal chromosomal locations. Complete probe sequences and the identity of the genes detected by the reference probes are available online (www.mrcholland.com).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mrcholland.com.

Length (nt)	Name
64-70-76-82 Q-fragments (only visible with <100 ng sample DNA)	
88-96	D-fragments (low signal indicates incomplete denaturation)
92	Benchmark fragment
100	X-fragment (X chromosome specific)
105	Y-fragment (Y chromosome specific)

MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mrcholland.com).

MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation ≤ 0.10 for all probes (except the *CEL* probes) over the experiment.

Required specimens

Extracted DNA, free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

Reference samples

A sufficient number (≥3) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from different unrelated individuals who are from families without a history of diabetes. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol (www.mrcholland.com).

Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (https://catalog.coriell.org) and Leibniz Institute DSMZ (https://www.dsmz.de/) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. Sample ID numbers from the Coriell Institute in the table below have been tested with this P357-A3 probemix at MRC Holland and can be used as positive control samples. The quality of cell lines can change; therefore samples should be validated before use.



Sample name	Altered target genes in P357-A3	Expected copy number alteration
NA00945	KLF11	Heterozygous deletion
NA01353	KLF11	Heterozygous duplication
NA04409	KLF11	Heterozygous duplication
NA10401	KLF11, NEUROD1	Heterozygous duplication
NA03918	NEUROD1	Heterozygous deletion
NA12519	PAX4	Homozygous duplication
NA13685	CEL	Heterozygous duplication
NA03330	PDX1	Heterozygous duplication
NA12606	PDX1	Heterozygous duplication
NA20359	HNF1B	Heterozygous duplication

Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mrcholland.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results

The standard deviation of each individual probe over all the reference samples should be ≤ 0.10 and the final ratio (FR) of each individual probe (except for the *CEL* probes) in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the FR of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

Copy number status	Final ratio (FR)
Normal	0.80 < FR < 1.20
Homozygous deletion	FR = 0
Heterozygous deletion	0.40 < FR < 0.65
Heterozygous duplication	1.30 < FR < 1.65
Heterozygous triplication/homozygous duplication	1.75 < FR < 2.15
Ambiguous copy number	All other values

Note: The term "dosage quotient", used in older product description versions, has been replaced by "final ratio" to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

- <u>Arranging probes</u> according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region or in or near the *KLF11* or *CEL* gene. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.

- <u>Normal copy number variation</u> in healthy individuals is described in the database of genomic variants: <u>http://dgv.tcag.ca/dgv/app/home</u>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- <u>Not all abnormalities detected by MLPA are pathogenic</u>. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- <u>Copy number changes detected by reference probes</u> or flanking probes are unlikely to have any relation to the condition tested for.
- False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more
 exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale
 peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net
 software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun
 the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount
 of sample by diluting PCR products.

P357 specific notes:

- *CEL* is located in a CNV region (http://dgv.tcag.ca/dgv/app/home). We observed duplications of this gene in several DNA samples from healthy individuals.
- *INS* gene mutations that cause MODY are all heterozygous missense mutations with a dominant-negative mode of action; a heterozygous deletion of the *INS* gene might not result in diabetes.

Limitations of the procedure

- In most populations, the major cause of genetic defects in MODY genes are small (point) mutations, most of which will not be detected by using SALSA MLPA Probemix P357 MODY Mix 2.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.

Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

KLF11, NEUROD1, PAX4, CEL, INS, PDX1 and HNF1B mutation databases

https://databases.lovd.nl/shared/genes/KLF11 https://databases.lovd.nl/shared/genes/NEUROD1 https://databases.lovd.nl/shared/genes/PAX4 https://databases.lovd.nl/shared/genes/CEL https://databases.lovd.nl/shared/genes/INS



https://databases.lovd.nl/shared/genes/PDX1 https://databases.lovd.nl/shared/genes/HNF1B

We strongly encourage users to deposit positive results in the Leiden Open Variation Database. Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on http://varnomen.hgvs.org/.

Please report copy number changes detected by the reference probes, false positive results due to SNVs and unusual results (e.g., a duplication of *INS* exons 1 and 3 but not exon 2) to MRC Holland: info@mrcholland.com.



Table 1. SALSA MLPA Probemix P357-A3 MODY Mix 2

Length		Chromosomal position (hg18) ^a					
(nt)	SALSA MLPA probe	Reference	HNF1B	PAX4/PDX1	KLF11/NEUROD1	CEL/INS	
64-105	Control fragments – see table in p	probemix cont	ent section f	or more information	ation	1	
131 *	Reference probe 16316-L22397	3q					
136 »	HNF1B probe 09858-L12290		Exon 7				
142	PAX4 probe 15491-L17346			Exon 1			
148 «	KLF11 probe 15492-L17347				Exon 2		
154	Reference probe 10195-L10655	16q					
160 »	HNF1B probe 07699-L12885		Exon 2				
165	PAX4 probe 15493-L17348			Exon 8			
178	Reference probe 05458-L04861	22q					
184	PAX4 probe 15495-L17350			Exon 3			
190 «	KLF11 probe 15496-L17351				Exon 4		
196	Reference probe 09102-L09161	4q					
202 »	HNF1B probe 08298-L18620		Exon 8				
208 »	INS probe 15497-L17352					Exon 2	
214 « ~	CEL probe 15498-L17353					Exon 2	
220	PDX1 probe 15499-L17354			Exon 2			
232 »	INS probe 15500-L17355					Exon 3	
238	HNF1B probe 07700-L07459		Exon 3				
244 »	INS probe 15501-L17356					Exon 1	
250 «	KLF11 probe 15502-L17357				Exon 3		
264 « ~	CEL probe 15503-L17358					Exon 7	
274	HNF1B probe 07702-L07461		Exon 4				
283	Reference probe 04404-L02610	14q					
292	PAX4 probe 15504-L17359			Exon 5			
301	HNF1B probe 08299-L09335		Exon 9				
309 « ~	CEL probe 15505-L17360					Exon 4	
317	NEUROD1 probe 15506-L17361				Exon 2		
326 « ~	CEL probe 15507-L18618					Exon 6	
335 »	HNF1B probe 07698-L18619		Exon 1				
346	Reference probe 10664-L11246	6р					
355	PDX1 probe 15508-L17363			Exon 1			
364 »	HNF1B probe 07701-L07460		Exon 4				
373	PAX4 probe 15509-L17364			Exon 6			
382	HNF1B probe 11551-L12298		Exon 5				
391	Reference probe 12522-L13572	18q					
400 « ~	CEL probe 15510-L17365			1	1	Exon 5	
409	PAX4 probe 15511-L17366			Exon 2	ĺ		
417	NEUROD1 probe 15512-L17367			1	Exon 1		
427	Reference probe 08787-L11328	10q		1	ĺ		
436	PAX4 probe 15513-L17368			Exon 7	1		
445 «	KLF11 probe 15514-L17369			1	Exon 1		
454	PAX4 probe 15515-L17370			Exon 4			
463 »	HNF1B probe 07704-L07463		Exon 6				
472	Reference probe 14846-L16554	3q					
481	PAX4 probe 15516-L17371	· · ·		Exon 9			
-	Reference probe 13814-L15308	5q				1	

^a See section Exon numbering on page 2 for more information.

* New in version A3.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

» The *INS* probes in this P357 probemix detect the same sequences as the *INS* probes in SALSA MLPA Probemix ME033. The indicated *HNF1B* probes in this probemix detect the same sequences as the *HNF1B* probes in SALSA MLPA Probemix P241.

~ *CEL* is located in a CNV region (http://dgv.tcag.ca/dgv/app/home). We observed duplications of this gene in several DNA samples from healthy individuals.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Table 2. P357-A3 probes arranged according to chromosomal location

Table 2a. KLF11 gene, 2p

Length (nt)	SALSA MLPA probe	KLF11 exon ^a	Ligation site NM_003597.5	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		start codon	168-170 (Exon 1)		
445 «	15514-L17369	Exon 1	183-184	ACACGCCGGACT-TCGCAGGCCCAG	2.5 kb
148 «	15492-L17347	Exon 2	286-287	TTGCAGCATCTT-GGAGCAGACAGA	2.0 kb
250 «	15502-L17357	Exon 3	1063-1064	AAGTAGCATGTT-ACCAGCTTTTTT	4.1 kb
190 «	15496-L17351	Exon 4	1553-1552 reverse	CTGCGCATGAAA-CGTCGGTCACAC	
		stop codon	1704-1706 (Exon 4)		

Table 2b. NEUROD1 gene, 2q

Length (nt)	SALSA MLPA probe	NEUROD1 exon ^a	Ligation site NM_002500.5	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		start codon	100-102 (Exon 2)		
417	15512-L17367	Exon 1	78-79	TCAGGACCTACT-AACAACAAAGGT	2.1 kb
317	15506-L17361	Exon 2	597-596 reverse	TGAACGAAGGAG-ACCAGGTCTGGG	
		stop codon	1168-1170 (Exon 2)		

Table 2c. PAX4 gene, 7q

Length (nt)	SALSA MLPA probe	PAX4 exon ^a	Ligation site NM_001366110.1	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		start codon	407-409 (Exon 3)		
142	15491-L17346	Exon 4 (1)	107 nt before exon 4 reverse	CCACCGAGTGCA-TCCTCTCCTGGG	0.6 kb
409	15511-L17366	Exon 5 (2)	606-607	CACAGGTGTCTT-GGAGCCAAAGGG	0.5 kb
184	15495-L17350	Exon 6 (3)	804-803 reverse	GTCCCTGGTCCT-CCTGTAATGCCC	0.6 kb
454	15515-L17370	Exon 7 (4)	30 nt before exon 7	TCACCATCACAA-ACCACCACAAAC	0.4 kb
292	15504-L17359	Exon 8 (5)	1 nt before exon 8	TCGATCTCCGCA-GAGTTCCAGCGT	0.5 kb
373	15509-L17364	Exon 9 (6)	1102-1101 reverse	TGCATTTCCCAC-TTGAGCTTCTCT	1.1 kb
436	15513-L17368	Exon 10 (7)	1168-1167 reverse	ACCTGTGCAGAG-ATGATTCCTGGG	0.4 kb
165	15493-L17348	Exon 11 (8)	1278-1279	AGAAAGGTGTCT-GAGTGACACCCC	0.6 kb
481	15516-L17371	Exon 12 (9)	1522-1521 reverse	GGCAGGACGGTA-AGGACAATGGGC	
		stop codon	1460-1462 (Exon 12)		

Table 2d. CEL gene, 9q

Length (nt)	SALSA MLPA probe	CEL exon ^a	Ligation site NM_001807.6	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		start codon	23-25 (Exon 1)		
214 « ~	15498-L17353	Exon 2	132-133	GGAAGGCGTCAA-TAAGAAGCTCGG	0.7 kb
309 « ~	15505-L17360	Exon 4	457-458	AACTACCTGTAT-GACGGCGAGGAG	1.5 kb
400 « ~	15510-L17365	Exon 5	645-644 reverse	CGAAGAGCGTGA-TGTTGTTGGGGT	0.3 kb
326 « ~	15507-L18618	Exon 6	732-733	CCGGCGAGCCAT-CAGCCAGAGCGG	0.2 kb
264 « ~	15503-L17358	Exon 7	823-824	GTGGGTTGCCCT-GTGGGTGATGCC	
		stop codon	2282-2284 (Exon 11)		



Table 2e. INS gene, 11p

Length (nt)	SALSA MLPA probe	INS exon ^a	Ligation site NM_000207.3	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		start codon	60-62 (Exon 2)		
244 »	15501-L17356	Exon 1	37-38	AAGAGGCCATCA-AGCAGGTCTGTT	0.4 kb
208 »	15497-L17352	Exon 2	233-232 reverse	TGCAGGTCCTCT-GCCTCCCGGCGG	1.0 kb
232 »	15500-L17355	Exon 3	10 nt after exon 3	GCCCTGCTGTGC-CGTCTGTGTGTC	
		stop codon	390-392 (Exon 3)		

Table 2f. PDX1 gene, 13q

Length (nt)	SALSA MLPA probe	PDX1 exon ^a	Ligation site NM_000209.4	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		start codon	140-142 (Exon 1)		
355	15508-L17363	Exon 1	190-191	TACAAGGACCCA-TGCGCGTTCCAG	5.7 kb
220	15499-L17354	Exon 2	2231-2232	TGTTCCGAGGTA-GAGGCTGTGCTG	
		stop codon	989-991 (Exon 2)		

Table 2g. HNF1B gene, 17q

Length (nt)	SALSA MLPA probe	HNF1B exon ^a	Ligation site NM_000458.4	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		start codon	176-178 (Exon 1)		
335 »	07698-L18619	Exon 1	180-181	TTGGAAAATGGT-GTCCAAGCTCAC	5.3 kb
160 »	07699-L12885	Exon 2	572-573	TGCAGCAACACA-ACATCCCCCAGA	5.8 kb
238	07700-L07459	Exon 3	768-769	AGACAAAAGCAG-TCAGGATCAGCT	2.0 kb
364 » #	07701-L07460	Exon 4	1051-1050 reverse	ACACGGACCTCA-GTGACCAAGTTG	0.3 kb
274	07702-L07461	Exon 4	73 nt after exon 4	CTTATCTGGTTT-AAGGGTTTTCAG	20.8 kb
382	11551-L12298	Exon 5	6 nt before exon 5	GTTTGTGTTGTT-TTGCAGGAGTGC	5.7 kb
463 »	07704-L07463	Exon 6	1399-1400	GTCTCAGGAGGA-GGTTTGCCCCCA	4.0 kb
136 »	09858-L12290	Exon 7	1634-1633 reverse	GCTCTGCTGCAT-GAGGGGCTGCTG	2.0 kb
202 »	08298-L18620	Exon 8	1798-1799	AGCAGCATCAGT-ACACTCACCAAC	12.0 kb
301	08299-L09335	Exon 9	2079-2080	CTCTCCCACGAT-GTCAAGGACTCC	
		stop codon	1847-1849 (Exon 9)		

^a See section Exon numbering on page 2 for more information.

^b Only partial probe sequences are shown. Complete probe sequences are available at www.mrcholland.com. Please notify us of any mistakes: info@mrcholland.com.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

» The *INS* probes in this P357 probemix detect the same sequences as the *INS* probes in SALSA MLPA Probemix ME033. The indicated *HNF1B* probes in this probemix detect the same sequences as the *HNF1B* probes in SALSA MLPA Probemix P241.

~ *CEL* is located in a CNV region (http://dgv.tcag.ca/dgv/app/home). We observed duplications of this gene in several DNA samples from healthy individuals.

This probe's specificity relies on a single nucleotide difference compared to a related gene or pseudogene. As a result, an apparent duplication of only this probe can be the result of a non-significant single nucleotide sequence change in the related gene or pseudogene.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Related SALSA MLPA probemixes

• ME033 TNDM: Contains probes for several genes involved in transient neonatal diabetes mellitus, including the *INS* probes that target the same sequences as the *INS* probes in this probemix (MODY 10) and *KCNJ11* probes (MODY 13).



- P117 ABCC8: Contains probes for the *ABCC8* gene, involved in familial hyperinsulinemic hypoglycemia 1 and linked to MODY 12.
- P241 MODY Mix 1: Contains probes for *HNF4A*, *GCK*, *HNF1A* and *HNF1B* genes, involved in MODY 1, 2, 3 and 5, respectively.

References

- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligationdependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.

Selected publications using SALSA MLPA Probemix P357 MODY Mix 2

- Horikawa Y et al. (2018). NEUROD1-deficient diabetes (MODY6): Identification of the first cases in Japanese and the clinical features. Pediatr Diabetes. 19:236-242.
- Morisada N et al. (2019). Clinical characteristics of HNF1B-related disorders in a Japanese population. Clin Exp Nephrol. 23:1119-1129.

P357 prod	P357 product history				
Version	Modification				
A3	One reference probe has been replaced.				
A2	Two reference probes and the 88 and 96 nt control fragments have been replaced (QDX2).				
A1	First release.				

Implemented changes in the product description

Version A3-01 – 28 July 2021 (04P)

- Product description rewritten and adapted to a new template.
- Ligation sites of the probes targeting all genes updated according to new versions of the NM_ reference sequences.
- Exon numbering of the PAX4 gene has been changed.
- Warning added to Table 2 for probe specificity relying on a single nucleotide difference between target gene and related gene or pseudogene.

Version 07 – 13 February 2018 (55)

- Product description adapted to a new product version (version number changed, lot number added, changes in Table 1 and Table 2, new picture included).
- ME033 added to related probemixes.
- Name changed from P357 MODY mix-2 to P357 MODY Mix 2.

Version 06 – 13 April 2017 (55)

- Warnings added in Table 1 and 2, 190 nt probe 15496-L17351, 214 nt probe 15498-L17353, 264 nt probe 15503-L17358, 309 nt probe 15505-L17360, 326 nt probe 15507-L18618, and 400 nt probe 15510-L17365.
- Minor textual changes.

More information: www.mrcholland.com; www.mrcholland.eu	
	MRC Holland bv; Willem Schoutenstraat 1
	1057 DL, Amsterdam, The Netherlands
E-mail	info@mrcholland.com (information & technical questions)
	order@mrcholland.com (orders)
Phone	+31 888 657 200