Product Description SALSA® MLPA® Probemix P405-B1 CMT1

To be used with the MLPA General Protocol.

Version B1

For complete product history see page 10.

Catalogue numbers:

- P405-025R: SALSA MLPA Probemix P405 CMT1, 25 reactions.
- P405-050R: SALSA MLPA Probemix P405 CMT1, 50 reactions.
- P405-100R: SALSA MLPA Probemix P405 CMT1, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see www.mrcholland.com).

Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mrcholland.com.

Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mrcholland.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

Intended purpose

The SALSA MLPA Probemix P405 CMT1 is an in vitro diagnostic (IVD)¹ or research use only (RUO) semi-quantitative assay² for the detection of deletions or duplications in the human *PMP22*, *MPZ* and *GJB1* genes in genomic DNA isolated from human peripheral whole blood specimens or buccal swabs. P405 CMT1 is intended to confirm a potential cause for and clinical diagnosis of Charcot-Marie-Tooth disease type 1 (CMT1), X-linked CMT (CMTX) or hereditary neuropathy with liability to pressure palsies (HNPP) and for molecular genetic testing of at-risk family members.

Copy number variations (CNVs) detected with P405 CMT1 should be confirmed with a different technique. In particular, CNVs detected by only a single probe always require confirmation by another method. Most defects in *MPZ* and *GJB1*, and some defects in *PMP22* are point mutations, none of which will be detected by MLPA. It is therefore recommended to use this assay in combination with sequence analysis.

Assay results are intended to be used in conjunction with other clinical and diagnostic findings, consistent with professional standards of practice, including confirmation by alternative methods, clinical genetic evaluation, and counselling, as appropriate. The results of this test should be interpreted by a clinical molecular geneticist or equivalent.

This device is not intended to be used for standalone diagnostic purposes, pre-implantation or prenatal testing, population screening, or for the detection of, or screening for, acquired or somatic genetic aberrations.

¹Please note that this probemix is for in vitro diagnostic (IVD) use in the countries specified at the end of this product description. In all other countries, the product is for research use only (RUO).

²To be used in combination with a SALSA MLPA Reagent Kit and Coffalyser.Net analysis software.

Clinical background

Charcot-Marie-Tooth disease (CMT), with a worldwide incidence of 1 in 2500, is the most common hereditary sensorimotor neuropathy, comprising a group of clinically and genetically heterogeneous peripheral neuropathies. CMT is characterized by progressive distal muscle atrophy and weakness, sensory disturbance, the absence of deep tendon reflexes, and *pes cavus* deformity of the foot. More than 80 different genes are



associated with CMT (http://www.ncbi.nlm.nih.gov/books/NBK1358/). Subtypes related to the genes *PMP22*, *GJB1*, *MPZ* and *MFN2* are the most common ones, being responsible for up to 95% of CMT cases with a final diagnosis (Padilha et al. 2020). The disease can be inherited in an autosomal dominant, autosomal recessive or X-linked manner. Table 1 provides an overview of the different genes involved in the CMT subtypes and the probemixes that cover these genes.

The most frequent form, CMT1A, accounts for as much as 50% of all CMT cases. CMT1A is a dominantly inherited, childhood-onset, slowly progressive motor and sensory neuropathy due to a duplication of *PMP22* on chromosome 17. CMT1B accounts for an additional ~10% of CMT1 patients. The gene involved in CMT1B is myelin protein zero (*MPZ*). CMT1A and 1B are clinically indistinguishable; classification is based solely on molecular findings.

X-linked CMT (CMTX) is characterized by moderate to severe motor and sensory neuropathy in affected males and usually mild to no symptoms in carrier females. Subtype CMTX1 accounts for about 90% of X-linked CMT, and is caused by mutations in the gap-junction protein beta 1 (*GJB1*) gene (https://www.ncbi.nlm.nih.gov/books/NBK1374/).

Hereditary neuropathy with liability to pressure palsies (HNPP) is characterized by repeated focal pressure neuropathies such as carpal tunnel syndrome and peroneal palsy with foot drop. Recovery from acute neuropathy is often complete; when recovery is not complete, the resulting disability is usually mild. Some affected individuals also have signs of a mild to moderate peripheral neuropathy. The prevalence of HNPP is estimated at 7-16 in 100,000. The penetrance is 100% but expressivity is highly variable even within the same family. Approximately 6-23% of individuals diagnosed with HNPP have an asymptomatic affected parent. A contiguous gene deletion of chromosome 17p12 that includes PMP22 is present in approximately 80% of affected individuals: 20% PMP22 the remaining have а pathogenic variant (https://www.ncbi.nlm.nih.gov/books/NBK1392/).

Table 1. Overview of the probemixes and genes related to CMT.

Probemix*	Genes and coverage	Condition	Remarks
P033-B4 CMT1 (IVD)	PMP22: all exons KIF1b: 2 probes	CMT1A and HNPP CMT2A1	PMP22 probes in P033-B4 have the same ligation site as <i>PMP22</i> probes in P405-B1 except for one exon 1 probe and one exon 4 probe. There is one additional PMP22 exon 5 probe present in P033-B4.
P405-B1 CMT1 (IVD)	PMP22: all exons MPZ: all exons GJB1: all exons	CMT1A and HNPP CMT1B CMTX	PMP22 probes in P405-B1 have the same ligation sites as PMP22 probes in P033-B4 except for one exon 1 probe and one exon 4 probe. There is one additional PMP22 exon 5 probe present in P033-B4. MPZ probes in P405-B1 have the same ligation sites as MPZ probes in P143-C3.
P143-C3 MFN2-MPZ (RUO)	MFN2: all exons MPZ: all exons	CMT2A CMT1B	MPZ probes in P143-C3 have the same ligation sites as MPZ probes in P405-B1.

^{*}IVD: in vitro diagnostic. RUO: research use only.

Gene structure

The *PMP22* gene spans ~35.5 kilobases (kb) on chromosome 17p12 and contains 5 exons. The *PMP22* LRG_263 is available at www.lrg-sequence.org and is identical to GenBank NG_007949.1.

The MPZ gene spans ~5.2 kb of genomic DNA on chromosome 1q23.3 and contains 6 exons. The MPZ LRG_256 is identical to GenBank NG_008055.1.





The *GJB1* gene spans ~10 kb of genomic DNA on chromosome Xq13.1 and contains 3 exons. The *GJB1* LRG_245 is identical to GenBank NG_008357.1.

Transcript variants

For *PMP22*, multiple variants have been described (https://www.ncbi.nlm.nih.gov/gene/5376). Transcript variant 1 encodes isoform 1 (NM_000304.4, 1828 nt, coding sequence 208-690). Alternative splicing results in multiple transcript variants. *PMP22* transcript variants 1-5 all encode isoform 1.

For MPZ, two variants have been described (https://www.ncbi.nlm.nih.gov/gene/4359). Transcript variant 1 encodes two isoforms, which result from the use of alternative in-frame translation termination codons. RefSeq transcript NM_000530.8 (1951 nt, coding sequence 64-810) represents the shorter isoform MPZ (also known as P0).

For *GJB1*, two variants have been described (https://www.ncbi.nlm.nih.gov/gene/2705). Transcript variant 1 represents the shorter transcript, and is transcribed from promoter P1 (NM_001097642.3, 1877 nt, coding sequence 69-920). Transcript variant 2 differs in the 5' UTR compared to variant 1, and is transcribed from promoter P2 (NM_000166.6; 1937 nt; coding sequence 129-980). Both variants 1 and 2 encode the same protein.

Exon numbering

The *PMP22* and *MPZ* exon numbering used in this P405-B1 CMT1 product description is the exon numbering from the LRG_263 and LRG_256 sequences, respectively. For *GJB1*, the exon numbering from the RefSeq transcript NM_001097642.3 is used. The exon numbering of the NM_ sequence that was used for determining a probe's ligation site does not always correspond to the exon numbering obtained from the LRG sequences. As changes to the databases can occur after release of this product description, the NM_ sequence and exon numbering may not be up-to-date.

Probemix content

The SALSA MLPA Probemix P405-B1 CMT1 contains 42 MLPA probes with amplification products between 130 and 445 nucleotides (nt). This includes 15 probes for the common 17p12 deletion/duplication region, two flanking probes for genes outside the common 17p12 deletion/duplication region, seven probes for the *MPZ* gene, and five probes for the *GJB1* gene. In addition, 10 ten reference probes are included that detect autosomal chromosomal locations and three probes detecting locations on the X-chromosome. Complete probe sequences and the identity of the genes detected by the reference probes are available online (www.mrcholland.com).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mrcholland.com.

Length (nt)	Name		
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)		
88-96	D-fragments (low signal indicates incomplete denaturation)		
92	Benchmark fragment		
100	X-fragment (X chromosome specific)		
105	Y-fragment (Y chromosome specific)		

MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mrcholland.com).





MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals of the same sex is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation ≤0.10 for all probes over the experiment.

Required specimens

Extracted DNA from human peripheral whole blood specimens or buccal swabs, free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

Reference samples

A sufficient number (≥3) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from different unrelated individuals who are from families without a history of CMT or HNPP. It is recommended to use samples of the same sex to facilitate interpretation. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol (www.mrcholland.com).

Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (https://catalog.coriell.org) and Leibniz Institute DSMZ (https://www.dsmz.de/) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. Sample ID numbers NA05167 and NA12214 from the Coriell Institute have been tested with this P405-B1 probemix at MRC Holland and can be used as a positive control samples to detect a heterozygous duplication of chromosome 17p12 that includes the *PMP22*, *COX10* and *TEKT3* genes. The quality of cell lines can change; therefore samples should be validated before use.

Performance characteristics

PMP22 duplications explain 50% of all CMT cases (https://www.ncbi.nlm.nih.gov/books/NBK1358/) and *PMP22* deletions explain 80% of all HNPP cases (https://www.ncbi.nlm.nih.gov/books/NBK1392/). *GJB1* deletions explain 0.3% of CMT cases (DiVincenzo et al. 2014). No duplications in *GJB1* have been described. No deletions or duplications in *MPZ* have been described in CMT cases, however, the association between *MPZ* mutations and CMT is well established (DiVincenzo et al. 2014). The analytical sensitivity and specificity for the detection of deletions or duplications in *PMP22*, *MPZ* and *GJB1* is very high and can be considered >99% (based on a 2006-2022 literature review).

Analytical performance can be compromised by: SNVs or other polymorphisms in the DNA target sequence, impurities in the DNA sample, incomplete DNA denaturation, the use of insufficient or too much sample DNA, the use of insufficient or unsuitable reference samples, problems with capillary electrophoresis or a poor data normalisation procedure and other technical errors. The MLPA General Protocol contains technical guidelines and information on data evaluation/normalisation.

Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mrcholland.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results

The expected results for the probes detecting autosomal sequences are allele copy numbers of 2 (normal), 1 (heterozygous deletion), or 3 (heterozygous duplication). The same results can be expected for the X-chromosome-specific probes in female samples. In rare cases, copy numbers of 0 (homozygous deletion) or





4 (heterozygous triplication/homozygous duplication) may be obtained. For the X-chromosome-specific probes in male samples, expected copy numbers are 1 (normal), 0 (deletion) or 2 (duplication).

The standard deviation of each individual probe over all the reference samples should be ≤ 0.10 and the final ratio (FR) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the FR of the probes can be used to interpret MLPA results when **reference samples of the same sex** have been used:

Copy number status		
Autosomal sequences and X chromosome sequences in females	X chromosome sequences in males	Final ratio (FR)
Normal	Normal	0.80 < FR < 1.20
Homozygous deletion	Deletion	FR = 0
Heterozygous deletion		0.40 < FR < 0.65
Heterozygous duplication		1.30 < FR < 1.65
Heterozygous triplication/homozygous duplication	Duplication	1.75 < FR < 2.15
Ambiguous copy number		All other values

Note: The term "dosage quotient", used in older product description versions, has been replaced by "final ratio" to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: http://dgv.tcag.ca/dgv/app/home. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for DMD by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- <u>Copy number changes detected by reference probes or flanking probes</u> are unlikely to have any relation to the condition tested for.
- False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.



P405 specific notes:

- Several types and subtypes of Charcot-Marie-Tooth disease exist, each involving a different gene or locus. Furthermore, there is a large overlap in clinical phenotype between the different types of CMT. Other disorders with comparable clinical features have been described, including hereditary neuralgic amyotrophy (HNA), amyloid neuropathies, Krabbe disease and hereditary ataxias. Therefore, molecular diagnosis may be complicated for certain patients.
- Please note that recurrent duplications have been described which can be detected by the two TEKT3 probes, but not by any of the PMP22 probes. These duplications may cause CMT through an unknown mechanism affecting PMP22 expression (Weterman et al. 2010; Zhang et al. 2010).

Limitations of the procedure

- In most populations, the major cause of genetic defects in the MPZ and GJB1 genes are small (point) mutations, none of which will be detected by using SALSA MLPA Probemix P405 CMT1.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect
 copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the
 possibility remains that biological changes in that gene or chromosomal region do exist but remain
 undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false
 positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe
 signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe
 oligonucleotide to the sample DNA.

Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

LOVD mutation database

https://databases.lovd.nl/shared/genes/PMP22 and https://databases.lovd.nl/shared/genes/MPZ and https://databases.lovd.nl/shared/genes/GJB1. We strongly encourage users to deposit positive results in the LOVD. Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on http://varnomen.hgvs.org/.

Please report copy number changes detected by the reference probes, false positive results due to SNVs and unusual results (e.g., a duplication of *PMP22* exons 2 and 4 but not exon 3) to MRC Holland: info@mrcholland.com.



Table 2. SALSA MLPA Probemix P405-B1 CMT1

Length (nt)	SALSA MLPA probe	Chromosomal position (hg18) ^a			
		Reference	MPZ	PMP22 / 17p12	GJB1
64-105	Control fragments – see table in pro	bemix content se	ection for more	information	
130	Reference probe 00797-L13645	5q			
136	GJB1 probe 02629-L02096				Exon 2
142	PMP22 probe 04656-L04039			Exon 1	
148	PMP22 probe 04657-L04461			Exon 3	
154	TEKT3 probe 22647-L32049			17p12, CMT1 region	
160	Reference probe 07394-L07041	12q			
165	MPZ probe 04895-L31858		Exon 1		
172	PMP22 probe 11539-L04463			Exon 4	
178	MPZ probe 04896-L04280		Exon 2		
184	PMP22 probe 02678-L02158			Upstream PMP22	
190	Reference probe 06148-L04992	7p			
196	GJB1 probe 13202-L14523				Exon 1
203	PMP22 probe 22645-L31860			Exon 3	
208	MPZ probe 04897-L32050		Exon 3		
213	GJB1 probe 13203-L19630				Exon 2
220	Reference probe 14968-L16704	6q			
226	MPZ probe 06139-L05583		Exon 1		
232	MPZ probe 04898-L17028		Exon 4		
238	PMP22 probe 22648-L31863			Exon 1	
244	PMP22 probe 04659-L19632			Exon 5	
250	Reference probe 18056-L31125	16q			
256	PMP22 probe 01462-L00927			Exon 2	
270	MPZ probe 04899-L20744		Exon 5		
278 ¬	DRC3 probe 01452-L20745			17p11, outside CMT1 region	
286	Reference probe 07737-L21372	20q			
293	TEKT3 probe 22649-L32051			17p12, CMT1 region	
301	MPZ probe 04900-L04284		Exon 6		
311	PMP22 probe 22650-L31865			Exon 4	
319	Reference probe 15385-L17792	3р			
328∫	X chromosome 13522-L14327				Xp11
346 ¬	ELAC2 probe 01466-L00917			17p12, outside CMT1 region	
355	PMP22 probe 22651-L31866			Upstream PMP22	
364	Reference probe 05953-L05397	2p			
373	PMP22 probe 02729-L02156			Upstream PMP22	
384∫	X chromosome 13750-L15237				Xp22
391	GJB1 probe 22646-L31861				Exon 2
400	COX10 probe 01468-L19633			17p12, CMT1 region	
409	Reference probe 02718-L00732	14q			
418∫	X chromosome 00820-L20737				Xq25
427	COX10 probe 01469-L20736			17p12, CMT1 region	
436	GJB1 probe 06188-L20735				Exon 2
445	Reference probe 12526-L13576	4q			

^a See section Exon numbering on page 3 for more information.

 $\int X\text{-chromosome}$ probe. Used for the determination of X-chromosome copy number.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

[¬] Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.





Table 3. P405-B1 probes arranged according to chromosomal location

Table 3a. MPZ gene

Length (nt)	SALSA MLPA probe	MPZ exon ^a	Ligation site NM_000530.8	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		stop codon	808-810 (Exon 6)		
301	04900-L04284	Exon 6	733-734	ATGCAATGCTGG-ACCACAGCAGAA	0.2 kb
270	04899-L20744	Exon 5	678-679	TTGCACAAGCCA-GGAAAGGACGCG	0.3 kb
232	04898-L17028	Exon 4	536-537	CGGGGTCGTTCT-GGGAGCTGTGAT	0.4 kb
208	04897-L32050	Exon 3	386-387	CCCTCGCTGGAA-GGATGGCTCCAT	0.5 kb
178	04896-L04280	Exon 2	217-218	TGCACTGCTCCT-TCTGGTCCAGTG	2.6 kb
165	04895-L31858	Exon 1	88-89	CTCCCTCATCCA-GCCCCAGCCCTA	0.2 kb
		start codon	64-66 (Exon 1)		
226	06139-L05583	Exon 1	85 nt before exon 1	CTGCACATGCCA-GGCTGCAATTGG	

Table 3b. PMP22 and 17p12 region

Length	SALSA MLPA	Gene/exon ^a	Ligation site	Partial sequence ^b (24 nt	Distance to
(nt)	probe	5 1.4.00	-	adjacent to ligation site)	next probe
346 ¬	01466-L00917	ELAC2		TTGGTCCTGAAT-GAGAACTGTGCC	1.2 M b
	r	Start of common	chromosome 17p12 del	etion/duplication region	1
427	01469-L20736	COX10	NM_001303.4; 1070-1071	CTCCTGGCAGTT-TCCTCATTTCAA	0.2 kb
400	01468-L19633	COX10	NM_001303.4; 1273-1274	CCATCAATGCGT-ACATCTCCTACC	1.0 M b
		PMP22	NM_000304.4		
		stop codon	688-690 (Exon 5)		
244	04659-L19632	Exon 5	672-673	ATCTATGTGATC-TTGCGGAAACGC	8.6 kb
311	22650-L31865	Exon 4	454-455	CTCTGTTCCTGT-TCTTCTGCCAAC	0.1 kb
172	11539-L04463	Exon 4	401-400 reverse	TGGTGGCCTGGA-CAGACTGCAGCC	19.5 kb
203	22645-L31860	Exon 3	367-368	TCCACCACTGTT-TCTCATCATCAC	0.1 kb
148	04657-L04461	Exon 3	310-309 reverse	GAGATCAGTTGC-GTGTCCATTGCC	1.5 kb
256	01462-L00927	Exon 2	229-230	TGTTGCTGAGTA-TCATCGTCCTCC	4.5 kb
		start codon	208-210 (Exon 2)		
238	22648-L31863	Exon 1	145-146	AGAAATCTGCTT-GGAAGAAGGGGT	0.1 kb
142	04656-L04039	Exon 1	25-26	ACCACCAGGGAA-CATCTCGGGGAG	2.3 kb
355	22651-L31866	Upstream	2.2 kb before exon 1	GGTGCTAGAAAT-AGCCAGTCTCAT	4.2 kb
373	02729-L02156	Upstream	6.4 kb before exon 1	GCCTCCATGGTT-AGAGACTAGAAT	5.9 kb
184	02678-L02158	Upstream	12.4 kb before exon 1	TGAAGAGCCCTT-GGATACGGAAGG	26.4 kb
154	22647-L32049	TEKT3	NM_031898.3; 1440-1441	ACACCATCCAGA-CCCTGCAGCAGC	27.4 kb
293	22649-L32051	TEKT3	NM_031898.3; 267-268	CCCACTCCAATT-TGACCCATAGCC	2.6 M b
		End of common	chromosome 17p12 dele	etion/duplication region	
278 ¬	01452-L20745	DRC3		CGGATCTCCAAG-ATCGACTCCCTG	
	•		•		•

Table 3c. GJB1

Length (nt)	SALSA MLPA probe	GJB1 exon ^a	Ligation site NM_001097642.3	<u>Partial</u> sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
196	13202-L14523	Exon 1	307 nt after exon 1 reverse	ACTTGGCCACCT-TGCACCTAGGGG	7.6 kb
391 +	22646-L31861	Exon 2	483 nt before exon 2	TAAAGCAGCATA-TGACTCCCCAGC	0.5 kb
		start codon	69-71 (Exon 2)		
136	02629-L02096	Exon 2	86-87	TGGACAGGTTTG-TACACCTTGCTC	0.4 kb
436	06188-L20735	Exon 2	447-448	ACAAGGTCCACA-TCTCAGGGACAC	0.7 kb
		stop codon	918-920 (Exon 2)		
213	13203-L19630	Exon 2	1129-1130	AGTGCTCAAGGT-TACTGGGAGTGT	

^a See section Exon numbering on page 3 for more information.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Related SALSA MLPA probemixes

P033 CMT1 Contains probes for the 17p12 chromosomal region and *KIF1b* gene.

P143 MFN2-MPZ Contains probes for the *MFN2* and *MPZ* genes.

P307 SEPT9 Contains probes for the SEPT9 gene involved in hereditary neuralgic amyotrophy.

P369 Smith-Magenis Contains probes for the 17p11.2 Smith-Magenis region.

References

- DiVincenzo C et al. (2014). The allelic spectrum of Charcot-Marie-Tooth disease in over 17,000 individuals with neuropathy. *Mol Genet Genomic Med*. 2(6):522-529.
- Padilha JPD et al. (2020). Diagnostic yield of targeted sequential and massive panel approaches for inherited neuropathies. Clin Genet. 98(2):185-190.
- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligationdependent probe amplification. Nucleic Acids Res. 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat*. 28:205.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. Anal Biochem. 421:799-801.
- Weterman MA et al. (2010). Copy number variation upstream of PMP22 in Charcot-Marie-Tooth disease. *Eur J Hum Genet*. 18(4):421-428.
- Zhang F et al. (2010). Mechanisms for nonrecurrent genomic rearrangements associated with CMT1A or HNPP: Rare CNVs as a cause for missing heritability. *Am J Hum Genet*. 86:892-903.

Selected publications using SALSA MLPA Probemix P405 CMT1

- Brown et al. (2014). Point mutation analysis of *PMP22* in patients referred for hereditary neuropathy with liability to pressure palsies. *Open J Genet.* 4(6):426-433.
- Gui et al. (2016). A new next-generation sequencing-based assay for concurrent preimplantation genetic diagnosis of Charcot-Marie-Tooth disease type 1A and aneuploidy screening. J Genet Genomics. 43(3):155-159.
- He et al. (2018). Clinical and genetic investigation in Chinese patients with demyelinating Charcot-Marie-Tooth disease. *J Peripher Nerv Syst.* 23(4):216-226.

^b Only partial probe sequences are shown. Complete probe sequences are available at www.mrcholland.com. Please notify us of any mistakes: info@mrcholland.com.

[¬] Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

⁺ The ligation site of this probe is located 15 nt before exon 1 of NM_000166.6.



- Kulshrestha et al. (2017). Deletion of P2 promoter of *GJB1* gene a cause of Charcot-Marie-Tooth disease. *Neuromuscul Disord*. 27:766-770.
- Lupo et al. (2015). Complexity of the hereditary motor and sensory neuropathies: Clinical and cellular characterization of the MPZ p.D90E mutation. J Child Neurol. 30(11):1544-1548.
- Milley et al. (2018). Genotypic and phenotypic spectrum of the most common causative genes of Charcot-Marie-Tooth disease in Hungarian patients. *Neuromuscul Disord*. 28(1):38-43.
- Nagappa et al. (2020). PMP22 gene-associated neuropathies: Phenotypic spectrum in a cohort from India. J Mol Neurosc. 70:778-789.
- Shahrizaila et al. (2014). X-linked Charcot-Marie-Tooth disease predominates in a cohort of multiethnic Malaysian patients. *Muscle Nerve*. 49(2):198-201.
- Wang et al. (2015). Clinical and genetic spectra in a series of Chinese patients with Charcot-Marie-Tooth disease. *Clin Chim Acta*. 451(pt B):263-270.

P405 produ	P405 product history		
Version	Modification		
B1	Four reference probes and three target probes have been replaced.		
A1	First release.		

Implemented changes in the product description

Version B1-03 – 27 January 2023 (04P)

- Table 1 adjusted; probemix versions added and remarks about PMP22 probes corrected.
- Number of exons in the GJB1 gene corrected in the Gene structure section.
- Information about transcript variants of *GJB1* updated in the Transcript variants section.
- Source for exon numbering of the *GJB1* gene corrected; exon numbering used is from RefSeq transcript NM_001097642.3 and not from the LRG_245 sequence.
- Chromosomal band of the 418 nt X chromosome probe 00820-L20737 corrected in Table 2.
- Ligation site of the 226 nt MPZ probe 06139-L05583 corrected in Table 3a.
- Remark added to Table 3 about the location of the ligation site of the 391 nt GJB1 probe 22646-L31861.
- Morocco removed as country with IVD status.
- Various minor textual and layout changes.

Version B1-02-06 May 2021 (04P)

- Product description rewritten and adapted to a new template.
- Various minor textual or layout changes.
- Ligation sites of the probes targeting the *GJB1* gene updated according to new version of the NM_ reference sequence.
- UK has been added to the list of countries in Europe that accept the CE mark.

Version B1-01 − 19 June 2020 (02P)

- Product description adapted to a new product version (version number changed, changes in Table 2 and Table 3).
- Product description rewritten and adapted to a new template.
- Colombia added as country with IVD status.
- Removed information about P129, this product is discontinued.
- New references added.

Version A1-03 - 25 April 2019 (04)

- Product description restructured and adapted to a new template.
- Various minor textual changes.
- Intended use updated.
- Table 1 added to provide an overview of the probemixes and genes related to CMT.
- Ligation sites of the probes targeting the MPZ, PMP22, TEKT3 and COX10 genes updated according to new version of the NM_ reference sequence.

Version A1-02 - 20 December 2018 (02)

- Regulatory status section updated to also include Morocco and Israel.





More infor	More information: www.mrcholland.com; www.mrcholland.eu			
***	MRC Holland bv; Willem Schoutenstraat 1			
	1057 DL, Amsterdam, The Netherlands			
E-mail	info@mrcholland.com (information & technical questions)			
	order@mrcholland.com (orders)			
Phone	+31 888 657 200			

IVD	EUROPE* CE ISRAEL COLOMBIA
RUO	ALL OTHER COUNTRIES

^{*}comprising EU (candidate) member states, members of the European Free Trade Association (EFTA) and the UK. The product is for RUO in all other European countries.