

Product Description SALSA® MLPA® Probemix P341-B4/P342-C1 PKHD1

To be used with the MLPA General Protocol.

P341 version B4. As compared to version B3, two reference probes have been replaced and one reference probe has been added. For complete product history see page 9.

P342 version C1. As compared to version B3, one probe for *PKHD1* has been removed and two reference probes have been replaced. For complete product history see page 9.

Catalogue numbers:

- **P341-025R:** SALSA MLPA Probemix P341 PKHD1 mix 1, 25 reactions.
- **P341-050R:** SALSA MLPA Probemix P341 PKHD1 mix 1, 50 reactions.
- **P341-100R:** SALSA MLPA Probemix P341 PKHD1 mix 1, 100 reactions.

- **P342-025R:** SALSA MLPA Probemix P342 PKHD1 mix 2, 25 reactions.
- **P342-050R:** SALSA MLPA Probemix P342 PKHD1 mix 2, 50 reactions.
- **P342-100R:** SALSA MLPA Probemix P342 PKHD1 mix 2, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see www.mlpa.com).

Certificate of Analysis: Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mlpa.com.

Precautions and warnings: For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mlpa.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

General information: The SALSA MLPA Probemix P341/P342 PKHD1 is a **research use only (RUO)** assay for the detection of deletions or duplications in the *PKHD1* gene, which is associated with autosomal recessive polycystic kidney disease (ARPKD).

ARPKD is a hereditary and severe form of polycystic kidney disease, affecting the kidneys and the hepatic biliary tract. The clinical spectrum is widely variable, with most cases presenting during infancy. The fetal phenotypic features classically include enlarged and echogenic kidneys, as well as oligohydramnios secondary to a poor urine output. Up to 50% of the affected neonates die shortly after birth, as a result of severe pulmonary hypoplasia and secondary respiratory insufficiency. In the subset that survives the perinatal period, morbidity and mortality are mainly related to severe systemic hypertension, renal insufficiency, and portal hypertension due to portal-tract fibrosis. Defects in the *PKHD1* gene are the cause of all typical forms of ARPKD.

More information is available at <https://www.ncbi.nlm.nih.gov/books/NBK1326/>.

This SALSA MLPA Probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.

Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene>

For NM_ mRNA reference sequences: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide>

Locus Reference Genomic (LRG) database: <http://www.lrg-sequence.org/>

Exon numbering: The *PKHD1* exon numbering used in this P341-B4/P342-C1 PKHD1 product description is the exon numbering from the RefSeq transcript NM_138694.4, which is identical to the NG_008753.1

sequence. The exon numbering and NM_ sequence used have been retrieved on 04/2020. As changes to the NCBI database can occur after release of this product description, exon numbering may not be up-to-date.

Probemix content: The SALSA MLPA Probemix P341-B4 PKHD1 mix 1 contains 44 MLPA probes with amplification products between 136 and 472 nucleotides (nt). This includes 35 probes for the *PKHD1* gene and nine reference probes that detect autosomal chromosomal locations. The SALSA MLPA Probemix P342-C1 PKHD1 mix 2 contains 44 MLPA probes with amplification products between 130 and 463 nucleotides (nt). This also includes 35 probes for the *PKHD1* gene and nine reference probes that detect autosomal chromosomal locations. Both probemixes together include 70 probes targeting 65 out of the 67 exons of the *PKHD1* gene, one probe for each exon of the gene with the exception of exon 17 and 44, and two probes are present for exons 27, 35, 36, 37, and 53. Complete probe sequences and the identity of the genes detected by the reference probes are available online (www.mlpa.com).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mlpa.com.

Length (nt)	Name
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)
88-96	D-fragments (low signal of 88 nt and 96 nt fragment indicates incomplete denaturation)
92	Benchmark fragment
100	X-fragment (X chromosome specific)
105	Y-fragment (Y chromosome specific)

MLPA technique: The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mlpa.com).

MLPA technique validation: Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation ≤ 0.10 for all probes over the experiment.

Required specimens: Extracted DNA free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

Reference samples: A sufficient number (≥ 3) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from unrelated individuals who are from families without a history of kidney diseases. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol.

Positive control DNA samples: MRC-Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (<https://catalog.coriell.org>) and Leibniz Institute DSMZ (<https://www.dsmz.de/home.html>) have a diverse collection of biological resources which may be used as a positive control DNA sample in your MLPA experiments. The quality of cell lines can change; therefore samples should be validated before use.

Data analysis: Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mlpa.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results: The standard deviation of each individual probe over all the reference samples should be ≤ 0.10 and the dosage quotient (DQ) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the DQ of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

Copy number status	Dosage quotient
Normal	$0.80 < DQ < 1.20$
Homozygous deletion	$DQ = 0$
Heterozygous deletion	$0.40 < DQ < 0.65$
Heterozygous duplication	$1.30 < DQ < 1.65$
Heterozygous triplication/Homozygous duplication	$1.75 < DQ < 2.15$
Ambiguous copy number	All other values

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Incomplete DNA denaturation (e.g. due to salt contamination) can lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes or flanking probes are unlikely to have any relation to the condition tested for.
- When running MLPA products, the capillary electrophoresis protocol may need optimization. False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: lower injection voltage / injection time settings, or a reduced amount of sample by diluting PCR products.

Limitations of the procedure:

- In most populations, the major cause of genetic defects in the *PKHD1* gene are small (point) mutations, most of which will not be detected by using SALSA MLPA Probemix P341/P342 PKHD1.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region do exist but remain undetected.
- Sequence changes (e.g. SNPs, point mutations, small indels) in the target sequence detected by a probe can cause false positive results. Mutations/SNPs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.

Confirmation of results: Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the

DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

PKHD1 mutation database: <https://databases.lovd.nl/shared/genes/PKHD1>. We strongly encourage users to deposit positive results in the Leiden Open Variation Database (LOVD). Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report copy number changes detected by the reference probes, false positive results due to SNPs and unusual results (e.g., a duplication of *PKHD1* exons 5 and 7 but not exon 6) to MRC-Holland: info@mlpa.com.

Table 1a. SALSA MLPA Probemix P341-B4 PKHD1 mix 1

Length (nt)	SALSA MLPA probe	Chromosomal position (hg18) ^a	
		Reference	PKHD1
64-105	Control fragments – see table in probemix content section for more information		
136	Reference probe 08030-L07811	11q24	
142	PKHD1 probe 10679-L11261		Exon 34
148	PKHD1 probe 10663-L11245		Exon 1
154	PKHD1 probe 10694-L11276		Exon 60
160	PKHD1 probe 20958-L11255		Exon 21
166	PKHD1 probe 10671-L11253		Exon 16
172	PKHD1 probe 12185-L13107		Exon 35
178	PKHD1 probe 10676-L11258		Exon 28
185	PKHD1 probe 10666-L11248		Exon 7
190	Reference probe 09953-L10412	17p13	
197	PKHD1 probe 10688-L11270		Exon 49
202	PKHD1 probe 10697-L12697		Exon 67
209	PKHD1 probe 10665-L11247		Exon 4
218	PKHD1 probe 10695-L24064		Exon 61
226	PKHD1 probe 10672-L11254		Exon 20
232	PKHD1 probe 20959-L29093		Exon 53
238	PKHD1 probe 10668-L11250		Exon 11
244	PKHD1 probe 10683-L11265		Exon 40
250 *	Reference probe 06387-L21633	8p23	
257	PKHD1 probe 10692-L11274		Exon 57
265	PKHD1 probe 10674-L11256		Exon 23
274	PKHD1 probe 10696-L11278		Exon 64
283	PKHD1 probe 10680-L11262		Exon 35
292 *	Reference probe 08722-L28962	9q21	
301	PKHD1 probe 10669-L11251		Exon 14
310	PKHD1 probe 10689-L11271		Exon 50
319	PKHD1 probe 10677-L11259		Exon 29
328	PKHD1 probe 10682-L11264		Exon 37
337	Reference probe 07722-L07432	7p13	
346	PKHD1 probe 10664-L11246		Exon 3
355	PKHD1 probe 10684-L11266		Exon 42
364	PKHD1 probe 10675-L11257		Exon 25
373	PKHD1 probe 10693-L11275		Exon 58
382	Reference probe 13055-L14238	15q14	
391 ±	PKHD1 probe 10691-L11273		Exon 55
400	PKHD1 probe 10670-L11252		Exon 15
409	PKHD1 probe 10681-L11263		Exon 36
418	PKHD1 probe 10687-L11269		Exon 46
427	Reference probe 10035-L11450	2q37	
436 +	PKHD1 probe 10678-L11260		Exon 33
445	PKHD1 probe 10667-L11249		Exon 8
454	PKHD1 probe 10685-L11267		Exon 43
463 *	Reference probe 18379-L23434	12p11	
472	Reference probe 12761-L13877	4q12	

a) See above section on exon numbering for more information.

* New in version B4.

± SNP rs150925674 could influence the probe signal. In case of apparent deletions, it is recommended to sequence the region targeted by this probe.

+ SNP rs770522674 could influence the probe signal. In case of apparent deletions, it is recommended to sequence the region targeted by this probe.

Table 1b. SALSA MLPA Probemix P342-C1 PKHD1 mix 2

Length (nt)	SALSA MLPA probe	Chromosomal position (hg18) ^a	
		Reference	PKHD1
64-105	Control fragments – see table in probemix content section for more information		
130	Reference probe 05169-L04550	8q24	
136	PKHD1 probe 10719-L11301		Exon 41
142	PKHD1 probe 10699-L11281		Exon 5
148	Reference probe 04730-L04147	7q21	
155	PKHD1 probe 10714-L11296		Exon 32
160	PKHD1 probe 10686-L11268		Exon 45
166	PKHD1 probe 10729-L24067		Exon 62
172	PKHD1 probe 10701-L11283		Exon 9
178	PKHD1 probe 10722-L11304		Exon 48
184	PKHD1 probe 10710-L11292		Exon 26
190 *	Reference probe 21401-L31194	3q22	
196	PKHD1 probe 10723-L11305		Exon 51
204	PKHD1 probe 10712-L29205		Exon 30
211	PKHD1 probe 10730-L29204		Exon 63
217	PKHD1 probe 10698-L29206		Exon 2
226	Reference probe 10433-L29207	9q34	
230	PKHD1 probe 10721-L29208		Exon 47
238	PKHD1 probe 10724-L11306		Exon 52
244	PKHD1 probe 10716-L11298		Exon 37
256	PKHD1 probe 10702-L11284		Exon 10
265	PKHD1 probe 10728-L11310		Exon 59
274	PKHD1 probe 10708-L11290		Exon 22
283	PKHD1 probe 10732-L11314		Exon 66
292	PKHD1 probe 11900-L12706		Exon 19
301	PKHD1 probe 10726-L11308		Exon 54
310	PKHD1 probe 10704-L11286		Exon 13
317	PKHD1 probe 11898-L24065		Exon 12
325	PKHD1 probe 10711-L24066		Exon 27
337	Reference probe 09937-L12248	8q13	
346	PKHD1 probe 10713-L11295		Exon 31
355	PKHD1 probe 10727-L11309		Exon 56
364 *	Reference probe 06348-L05863	1p21	
373	PKHD1 probe 10718-L11300		Exon 39
382	PKHD1 probe 11901-L12707		Exon 27
391	PKHD1 probe 10725-L11307		Exon 53
400	PKHD1 probe 20685-L22083		Exon 36
409	Reference probe 09999-L10331	20q13	
416	PKHD1 probe 20960-L29094		Exon 6
422	PKHD1 probe 12187-L29298		Exon 18
427	Reference probe 15893-L24326	2p16	
437	PKHD1 probe 10731-L11313		Exon 65
445	PKHD1 probe 10709-L11291		Exon 24
454	PKHD1 probe 10717-L11299		Exon 38
463	Reference probe 09908-L10321	16p13	

a) See above section on exon numbering for more information.

* New in version C1.

Table 2. PKHD1 probes arranged according to chromosomal location

Length (nt)	SALSA MLPA probe	PKHD1 exon ^a	Ligation site NM_138694.4	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		<i>start codon</i>	<i>265-267 (Exon 2)</i>		
148	10663-L11245	Exon 1	81-82	CTCTAACCCAGAT-AACATGTCCACG	2.6 kb
217	10698-L29206	Exon 2	296-297	TCTGATGAGTAT-TGAAGTACTACT	1.7 kb
346	10664-L11246	Exon 3	348-349	ATTGAACCTGAA-GAAGGTAGCCTT	0.7 kb
209	10665-L11247	Exon 4	436-437	ATGGCTCTCAAT-TGGAGATACACC	2.5 kb
142	10699-L11281	Exon 5	580-581	AGGGTCTGTACT-TCCTGGAAGCAT	3.7 kb
416	20960-L29094	Exon 6	700-701	TTTATCCACCAA-GTGGTGTCCAG	2.8 kb
185	10666-L11248	Exon 7	759-760	GGAAGATTGGAA-ACTTTTGATTTT	1.3 kb
445	10667-L11249	Exon 8	821-822	TCAAGGAGACAA-ATGGGTTACTCC	1.1 kb
172	10701-L11283	Exon 9	8 nt before exon 9	TCAGTTTGTCTT-TTCCTTAGTTAT	0.7 kb
256	10702-L11284	Exon 10	959-960	CTTCTCAGTATT-TAACAAAGGAAA	0.9 kb
238	10668-L11250	Exon 11	1001-1002	GGCATGGCTGAT-CAGTGCTAAACA	3.4 kb
317	11898-L24065	Exon 12	1065-1066	GTGTTTCCAGAA-ACTGGGAGCCTT	1.0 kb
310	10704-L11286	Exon 13	1184-1185	TCCCAGGAAGAT-TGAGTGCACCAC	2.5 kb
301	10669-L11251	Exon 14	1349-1350	CAGTTCCTCATT-TGGGTTTTGGTC	2.6 kb
400	10670-L11252	Exon 15	1479-1480	TTCAGTTGGTCA-GAGGAACCAAGG	1.5 kb
166	10671-L11253	Exon 16	1699-1700	ACACCTGGCTGA-ATCCTGATGTGG	1.7 kb
	No probe	Exon 17			
422	12187-L29298	Exon 18	1906-1907	CAGTAAAATGCA-AACTGGAACCCC	1.2 kb
292	11900-L12706	Exon 19	61 nt after exon 19	TGCTCTGGAT-TCAAGACTGAAA	1.4 kb
226	10672-L11254	Exon 20	2128-2129	AAGCCACATGA-ACAAGATCCTGA	0.9 kb
160	20958-L11255	Exon 21	2263-2264	TCTGGGAGACTT-GTGTGCGTTGCT	2.9 kb
274	10708-L11290	Exon 22	2418-2419	TCTCAAGCTGAT-TCTGGAACGGCT	1.7 kb
265	10674-L11256	Exon 23	2558-2559	TGTGCCCACTGA-AGGAACAGAAGA	2.4 kb
445	10709-L11291	Exon 24	2688-2689	CCTGTACAAATT-TCTGCTCATCAC	1.2 kb
364	10675-L11257	Exon 25	2926-2927	ATGGTGGAGTTT-TTCTTGACCCA	1.3 kb
184	10710-L11292	Exon 26	3030-3031	GGTTCCTGCTCT-TTCCAGTACCTC	0.6 kb
382	11901-L12707	Exon 27	3167-3168	TACAGTGAACAA-AACGAGTTGCAA	0.1 kb
325	10711-L24066	Exon 27	3254-3255	TCGGATCTTGAT-GTTGGTGAGACC	7.4 kb
178	10676-L11258	Exon 28	3459-3460	GCTACAAGCAAT-TCAAGCAGAATT	2.5 kb
319	10677-L11259	Exon 29	3522-3523	GTGAATGTGACT-GTGATCAGAGGG	4.8 kb
204	10712-L29205	Exon 30	3679-3680	ACTATACGGATT-TGGATGTGGAAG	0.4 kb
346	10713-L11295	Exon 31	3857-3858	CCTCACAGAAGT-TTTCAGCATCGA	2.4 kb
155	10714-L11296	Exon 32	4562-4563	TACTTGTGTGAT-TTTGAGTTGGG	2.7 kb
436 +	10678-L11260	Exon 33	5604-5605	CTGGCTAATGCT-ACAGTGTCTGCC	5.4 kb
142	10679-L11261	Exon 34	5790-5791	ATTTGCGAGGAA-AGTTCCAATGC	7.0 kb
172	12185-L13107	Exon 35	22 nt before exon 35	ATGACACCCCAT-TTAACCTCCCCT	0.1 kb
283	10680-L11262	Exon 35	5966-5967	CAATCAGCCAAT-TACCGTCAAGAT	50.4 kb
409	10681-L11263	Exon 36	6068-6069	TCACAGCTGGTT-TCCTGAAAGGCT	0.1 kb
400	20685-L22083	Exon 36	6126-6127	GCCCAATTGCTT-CTGCTGGCACT	25.7 kb
244	10716-L11298	Exon 37	6234-6235	GCCATCCTTGTT-TCTGATGGTGGA	0.1 kb
328	10682-L11264	Exon 37	6302-6303	TCAGATCACACT-CTACGGGAGTTC	21.6 kb
454	10717-L11299	Exon 38	6423-6424	TGTCTTAGAGCA-ACTGCCCATGCC	0.7 kb
373	10718-L11300	Exon 39	6680-6681	CAGCAGGAGTAT-TACCATAACAAGG	2.5 kb
244	10683-L11265	Exon 40	6882-6883	GGAGTGCAGTTT-CAAGTCTTGGGG	3.0 kb
136	10719-L11301	Exon 41	6969-6970	CAGGGCTGCACA-GTGAGGAACTCC	2.3 kb
355	10684-L11266	Exon 42	7095-7096	GAGATGAGATAT-ATCTCCTGGGAG	0.4 kb
454	10685-L11267	Exon 43	7233-7234	TCTGGCATCTAT-ATCTGCAGTCCC	17.7 kb
	No probe	Exon 44			
160	10686-L11268	Exon 45	7448-7449	GTTCCAGAGCTT-CACAGTTTGGGA	2.8 kb
418	10687-L11269	Exon 46	7577-7578	AAATACTTCAGT-TACTGACAGCTT	12.5 kb
230	10721-L29208	Exon 47	7643-7644	ATCTGGGATTAA-AACTCCTAAAAG	2.5 kb
178	10722-L11304	Exon 48	7789-7790	TGAAGTTTACAA-ACTCTTCAAACCT	12.1 kb
197	10688-L11270	Exon 49	8124-8125	CTCTTGGACCAA-GAGACCTACTCA	8.1 kb
310	10689-L11271	Exon 50	8296-8297	TGTCTTTTCAT-TTCTTCCATCAC	11.4 kb
196	10723-L11305	Exon 51	8399-8400	AGTTCAAGTCAT-TCTCCGGGTGAA	5.5 kb
238	10724-L11306	Exon 52	8526-8527	GGATACAACAAT-ACCATTCCAGGC	39.6 kb

Length (nt)	SALSA MLPA probe	PKHD1 exon ^a	Ligation site NM_138694.4	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
232	20959-L29093	Exon 53	8594-8595	GGATACAGATCT-TCCATTCTTCAA	0.1 kb
391	10725-L11307	Exon 53	8675-8676	TGTGGCATGCAT-GGTCATTGCAGG	15.4 kb
301	10726-L11308	Exon 54	8774-8775	AGAGGGAGTCTT-TTGTGACCGTAT	3.1 kb
391 ±	10691-L11273	Exon 55	8845-8846	TTCATCTTTACA-GTGCTTATCCTA	17.9 kb
355	10727-L11309	Exon 56	8985-8986	GAGCCTCATGAA-GCAGAGGTCTC	1.6 kb
257	10692-L11274	Exon 57	9134-9135	GACCCGAAATAT-ACAAATTCAGCC	4.7 kb
373	10693-L11275	Exon 58	9273-9274	TTGTAATCATCT-GTTGAATTCAGT	1.8 kb
265	10728-L11310	Exon 59	10161-10162	ATTCTACCAAAT-GCAGAGAACAGT	2.4 kb
154	10694-L11276	Exon 60	10375-10376	CAGTTTCTGTAT-TTCCTAAAACAG	84.9 kb
218	10695-L24064	Exon 61	10871-10872	TGCCAACTATTT-CAACATCATGGA	10.4 kb
166	10729-L24067	Exon 62	11516-11517	CCAGCCTTCAGA-TGGAGAAGTGGG	1.1 kb
211	10730-L29204	Exon 63	11637-11638	GCTTCCCTGGAA-GGAGCATCAGAC	9.1 kb
274	10696-L11278	Exon 64	11685-11686	GCAGAACTCAA-GATGGTTATGTT	6.3 kb
437	10731-L11313	Exon 65	11825-11826	TGTGACTAGGAA-GGAGAAGTCGAC	5.7 kb
283	10732-L11314	Exon 66	7 nt after exon 66	AGAAGGTAAGCT-TGGAGGGTGGAG	7.7 kb
202	10697-L12697	Exon 67	12307-12308	TGGCAGGCCAAA-ATCAGCTGCTGC	
		stop codon	12487-12489 (Exon 67)		

a) See above section on exon numbering for more information.

b) Only partial probe sequences are shown. Complete probe sequences are available at www.mlpa.com. Please notify us of any mistakes: info@mlpa.com.

± SNP rs150925674 could influence the probe signal. In case of apparent deletions, it is recommended to sequence the region targeted by this probe.

+ SNP rs770522674 could influence the probe signal. In case of apparent deletions, it is recommended to sequence the region targeted by this probe.

References

- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.


Selected publications using SALSA MLPA Probemix P341/P342 PKHD1

- Koczok K et al. (2018). Interfering effect of maternal cell contamination on invasive prenatal molecular genetic testing. *Prenatal diagnosis*, 38(9), 713-719.
- Lazaros L et al. (2019). Identification of a Novel Intragenic Deletion of the PKHD1 Gene in a Patient with Autosomal Recessive Polycystic Kidney Disease. *International Journal of Pediatrics*, 7(10), 10291-10297.
- Melchionda S et al. (2016). Expanding the mutation spectrum in 130 probands with ARPKD: identification of 62 novel PKHD1 mutations by sanger sequencing and MLPA analysis. *J Hum Genet.* Epub 2016 May 26. doi: 10.1038/jhg.2016.58.
- Miyazaki J et al. (2015). Intragenic duplication in the PKHD1 gene in autosomal recessive polycystic kidney disease. *BMC Med Genet.* 16(98): 2-6.
- Obeidova L et al. (2015). Molecular genetic analysis of PKHD1 by next-generation sequencing in Czech families with autosomal recessive polycystic kidney disease. *BMC Med Genet.* 16(116):1-12.
- Szabó T et al. (2018). Comprehensive genetic testing in children with a clinical diagnosis of ARPKD identifies phenocopies. *Pediatric Nephrology*, 33(10), 1713-1721.
- Zvereff V et al. (2010). Identification of PKHD1 multiexon deletions using multiplex ligation-dependent probe amplification and quantitative polymerase chain reaction. *Genet Test Mol Biomarkers.* 14:505-10.

P341 Product history	
<i>Version</i>	<i>Modification</i>
B4	Two reference probes have been replaced and one reference probe has been added.
B3	The lengths of several probes have been adjusted.
B2	Two reference probes have been removed and control fragments have been adjusted (QDX2).
B1	Extra exon 35 probe, X and Y fragments added, and several reference probes have been replaced.
A1	First release.

P342 Product history	
<i>Version</i>	<i>Modification</i>
C1	One probe for <i>PKHD1</i> has been removed and two reference probes have been replaced.
B3	The lengths of several probes have been adjusted.
B2	One reference probe has been replaced, one additional reference probe has been included, and the control fragments have been adjusted (QDX2).
B1	Exon 18 probe added, exon 17 probe removed, X and Y fragments added, and one reference probe removed.
A1	First release.

Implemented changes in the product description	
<i>Version B4/C1-01 — 12 June 2020 (02P)</i>	
<ul style="list-style-type: none"> - Product description rewritten and adapted to a new template. - Product description adapted to a new product version (version number changed, changes in Table 1 and Table 2). - Ligation sites of the probes targeting the <i>PKHD1</i> gene updated according to new version of the NM_ reference sequence. - Small changes of probe lengths in Table 1 and 2 in order to better reflect the true lengths of the amplification products. - Warning added to Tables 1 and 2 about specific SNPs that could influence the probe signal. 	
<i>Version 08 (55)- 28 June 2016</i>	
<ul style="list-style-type: none"> - Product description adapted to a new product version (version number changed, lot number added, small changes in Table 1 and Table 2, new picture included). - Various minor textual changes. - New references added on page 2. - <i>PKHD1</i> exon numbering adjusted. 	
<i>Version 07 (49)</i>	
<ul style="list-style-type: none"> - Figure(s) based on the use of old MLPA buffer (replaced in December 2012) removed. 	
<i>Version 06 (49)</i>	
<ul style="list-style-type: none"> - Product description adapted to a new lot (lot number added, small changes in Table 1 and Table 2, new picture included). 	
<i>Version 05 (48)</i>	
<ul style="list-style-type: none"> - Electropherogram pictures using the new MLPA buffer (introduced in December 2012) added. 	

More information: www.mlpa.com; www.mlpa.eu	
	MRC-Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands
E-mail	info@mlpa.com (information & technical questions); order@mlpa.com (orders)
Phone	+31 888 657 200